

# NEWA PLUS WIRELESS EDITION



Quick Guide

## Before you start

- Make sure your NEWA+ is fully charged before starting treatments. We recommend charging NEWA+ overnight before starting a treatment.
- Remove the charger from the NEWA+ device
- You can treat your forearm first to see how your skin reacts to NEWA+.
- Remove all jewelry from the treatment areas. Thoroughly clean and dry your face.
- Choose your treatment area (from a total of 3 areas).
- It is recommended to perform 3-5 treatments a week for a period of one month. For maintenance use NEWA+ twice a week.



## Start your NEWA+ treatment:



### Step 1

Remove the cover and apply two thick lines of NEWA gel, making sure the electrodes area is covered completely.



### Step 2

Spread the NEWA gel on the treatment area in a circular motion.



### Step 3

Turn on the device by moving the power level switch to 2. The blue LED will turn on.



### Step 4

Press and release the button. The blue LED will start blinking.



### Step 5

Begin treatment. Make sure the active tip area and electrodes are in full contact with the skin. The LED will blink green when your skin reaches optimal treatment temperature. If you find treatment at level 2 uncomfortable, move the power switch to level 1.

## Treatment

- Do not stop the circular movement of the device.
- NEWA+ will vibrate after 4 minutes signaling the end of the treatment. Apply more gel before moving to the next treatment area.
- At the end of the treatment gently clean your face and remove excess gel.
- If NEWA+ is turned on but not in use, it will automatically turn off after a few minutes. To restart the NEWA+, move the power switch to 0, and then slide it back to 2.
- To maximize battery life, it is recommended to charge NEWA+ at least once every two weeks, over night.
- To get the best out of your NEWA device, it is vital to properly clean it after each treatment. Simply follow these easy instructions:
  1. Make sure your NEWA is switched off and disconnected from the power source.
  2. Lightly dampen a washcloth with water.
  3. Gently wipe the electrodes of your device to remove any remaining gel residue.
  4. Gently wipe down the case of your device.
  5. Let your NEWA air dry

## NEWA+ LED indications:

LED Color	Description
 Solid Orange and Green LED	Device is plugged in and is charging (device is off).
 Solid Blue LED	Device is on and ready for treatment
 Blinking Blue LED	Device is delivering RF energy, and has not yet reached the optimal treatment temperature.
 Blinking Green LED	Indicates that the optimal treatment temperature has been achieved.
 Blinking Orange LED	Battery is running low during treatment. Continue with your full treatment session. When the treatment is complete, connect the device to the charger.
 Solid Green LED	Device is fully charged (plugged in) and is turned off.
 No LED Indication	Device is turned off and not plugged in.

For more information check the UM or visit [www.mynewa.com](http://www.mynewa.com)

NEWA  PLUS  
WIRELESS  
EDITION